



One ml contains 100 mg  
100 mg/ml  
Testosterone Base

**TEST BASE<sup>®</sup>**



10 ml vial



**TEST BASE<sup>®</sup>**

Testosterone Base 100 mg/ml

**COMPOSITION:**

Each ml of the solution contains:  
Testosterone Base - 100 mg,  
Benzyl Alcohol, see the package  
for further information.  
Solution for injection. Read the  
package Leaflet before use.

**Manufacturer:**

Canadian Pharmaceuticals  
32 Tower Ave  
Steinbach, MB R5G 0B8  
Canada  
[www.thecanadianpharm.com](http://www.thecanadianpharm.com)

10 ml vial



**TEST BASE<sup>®</sup>**

Testosterone Base 100 mg/ml

**How to store TEST BASE**

It should be stored in a dry place away  
from heat and humidity, 15° to 25° C (59°  
to 77° F).

Keep Test Base out of the reach of  
children and away from pets.

Store in the original package in order to  
protect from light. Do not throw away any  
medicines via wastewater or household  
waste. Ask your Pharmacist how to throw  
away medicines you no longer use.  
MEDICAL PRODUCT SUBJECT TO  
MEDICAL PRESCRIPTION.



LOT L80546738  
EXP 10.2023